




**DIA-STRON**  
DELIVERING MEASUREMENT SOLUTIONS

# CYC802 Cyclic Module



## Overview

The CYC802 Cyclic Tester simulates a realistic approach to everyday hair grooming by subjecting fibres to repeated cyclic tensile deformations until failure. Data analysis from the CYC802 shows greater differences between samples than conventional tensile stress/strain experiments.

Principal benefits:

- High throughput testing
- High capacity: 50 sample linear cassette and up to 4 cassettes in a cassette hotel
- Automated operations and analysis
- Up to 4 cyclic modules can be integrated on an automated platform
- Multi-tasking allows for simultaneous measurement, reducing resting time

Applications and claims:

- Fatigue testing using repeated strain, force or stress
- Hair root vs tip survival analysis
- Impact of hair treatments
- S-N curves



### **Metrology principle —**

The CYC802 fatigue tester is designed around a voice-coil drive which repeatedly brings the sample to a pre-determined strain, force or stress. The sample is mounted using Dia-Stron brass crimps and placed within two sample pockets. With the fibre dimensions already captured, the CYC802 brings the sample to a pre-set force before starting the fatigue measurement. Various factors contributing to fibre failure are: presence and propagation rate of flaws depending on ethnicity, chemical or physical damage, grooming regime and environmental factors such as UV exposure.

### **Dedicated software – UvWin**

The CYC802 system is controlled using Dia-Stron UvWin software. Parameters for these methods can be easily edited within the software. UvWin enables automatic data correction for system compliance.

UvWin also offers a number of integrated data processing tools and the raw data can be exported as a text file, for use in Excel or other statistical packages. UvWin is compatible with the latest versions of the Windows OS.

# Sample data and analysis —

Figure 1: First cycle









Figure 2: Failure cycle



Kaplan-Meier plot using a linear x-axis (left) and a logarithmic x-axis (right)



Boxplots of elastic moduli  
Blue dots represent means



Scatterplots of elastic moduli  
Blue dots represent means  
Red dots represent medians



## References —

### Publications:

Marsh JM., Mamak M., Wireko F., Lebron A., Cambron T., Huber T., David IN., Williams REA. and McComb DW. (2018): “Multimodal Evidence of Mesostructured Calcium Fatty Acid Deposits in Human Hair & Their Role on Hair Properties”; ResearchGate, September 2019

Everaert EPJM., Zhang S., Tran D., Kroon B., Zhang G., Thompson B. and McMullen RL. (2015): “Strengthening the Hair Fiber from Within: Repairing the Cortex of Damaged Hair”; IFSCC Conference, Zurich 2015

Evans TA. (2009): “Fatigue testing of hair—A statistical approach to hair breakage”; Journal of Cosmetic Science, 2009, 60, 599-616

### Examples of use in patent claims:

US20170007518A1 Method of strengthening hair fibers and protecting dyed hair color from fading or wash-out (Ashland) May 2019

### Uses for claims in technology advertising:

- Ashland
- Evonik

### CYC802 Module

Extension range	28-50mm
Speed range	1-100mm/sec
Force range	0 to 20N (2000gmf)
Force resolution	0.05gmf
Displacement resolution	10 $\mu$ m
Dimensions (LxWxH)	445x20x130mm
Weight	2.25kg

### Programmable Features

Methods	Constant force Constant strain Constant stress S-N curve
---------	---

### Content

UV1000  
PU1100  
CYC802 module  
USB and power cables  
UvWin software for Windows OS

### Requirements

Power supply	Universal 85-265V AC 47-63Hz, 100W
Computer	Windows OS: 7 and 10 1 x USB port

# Contact Us

## **Dia-Stron Ltd.**

9 Focus Way  
Andover, Hampshire  
SP10 5NY | United Kingdom  
T. +44(0) 1264 334700

## **Dia-Stron Inc.**

9 Trenton Lakewood Road  
Clarksburg, NJ  
08510 | U.S.A.  
T. +1 (609) 454 6008

Email: [enquiry@diastron.com](mailto:enquiry@diastron.com)

[www.diastron.com](http://www.diastron.com)